

# WORLD-CLASS EDUCATION

The range of applications using diode laser continue to grow, in all fields of medicine.

Through science and operational-based educational training, laser treatment can be easily applied.

Medency provides laser certification with the support of the best world-class Academies and Universities in the field.

Increase the confidence with new technologies, establish your modern facilities.

## RAPIDO Surgical Laser System

Code  
Wavelength / power  
Frequency  
Pulse mode  
Pulse range  
Medical /Laser class  
Dimensions/ Weight

**LADER102.2**  
**980nm @10W**  
**Up to 25 kHz**  
**Continuous - Multiple pulse**  
**20µs up to 1s**  
**IIB / 4**  
**205x130x50mm**

Rapido is a CE medical device with CE trademark. In accordance with EEC Directives for Medical Device: EEC Directives 93/42. All Medency's products are made in Italy.



**MEDENCY**  
TECHNOLOGY | INNOVATION | PASSION

# THE MULTIDISCIPLINARY SURGICAL LASER




**RAPIDO  
MED**



**MEDENCY Srl Piazza della Libertá 49**  
**I-36077 Altavilla - Vicenza Italy**  
**Tel +39 0444 371462**

[info@medency.com](mailto:info@medency.com) | [www.medency.com](http://www.medency.com)



YOUR LOCAL DEALER

For medical professionals only ©2020 MEDENCY. All rights reserved. Code DPSURE01

[medency.com](http://medency.com)



Diode Laser Technology for advanced medical procedures

# SAFE LASER ENERGY FOR MINIMAL INVASIVE SURGERY

What makes Rapido Med diode laser so versatile is its 980 nm wavelength, which allows its use in almost any medical surgical procedure, such as:

#### PLDD

The Percutaneous Laser Disc Decompression is a minimally invasive technique, which safe, time-effective and can be performed in a total out-patient setting.

#### EVLV

The action of laser light offers an efficient and non-invasive action for treating the source of varicose veins, providing patients with least recovery time and a fast return to normal daily routine.

#### GYNAECOLOGY

The therapeutic benefits of the laser light applied on the tissue bring a tightening effect of vaginal wall and allow the reconstruction of the weakened tissues responsible for vaginal relaxation (LVR technique).

#### PROCTOLOGY

The action of diode laser provides fast cutting with immediate coagulation, limiting damages on surrounding tissues. The main applications which find reliable and effective results are hemorrhoidectomy and rectal fistula removal.



**RAPIDO**  
MED

Medency is constantly up-to-date on the latest clinical findings and focused on developing breakthrough laser technology to increase the reliability of the studio and patients' experience.

## YOUR PATIENTS DESERVE THE BEST!

- Effective surgery
- Minimal thermal effect
- High working efficiency
- Faster and less invasive treatments
- Reduced recovery time and costs
- No post OP side effects
- Easy-to-use software
- Portable unit working with battery and wireless pedal
- Ergonomic and attractive design
- Flexible fiber and versatile accessories

## WIDE RANGE OF OPTICAL FIBERS



Diabetic Foot



PLDD



EVLV



DCR Surgery



Gynaecology



Proctology